

KN 95

ANTI PARTICULATE FACE MASK

Model: AC9503

GB2626-2006



Refreshing
and not stuffy



Fog prevention
and haze



Quiet enjoyment



Small and light



50 pcs / box
2 pcs / bag

PACKAGING

Product



2 pcs/Bag




25bags(50pcs)/Box



TECHNICAL DATASHEET

1. Anti Particulate Face Mask, 4-ply



2. Size: (see photos)



3. Material:

- (1) Non-woven Material 3 pcs
- (2) Meltblown Material 1 pcs
- (3) Aluminum nosepiece 1 pcs
- (4) Nylon & Spandex Webbing 2 pcs

4. Package:

50pcs/box, 20boxes/ctn

CE & TEST REPORT

Certificate - Сертификат - 證明書 - Certificat - 증명서 - شهادة



Certificate of Compliance

No. CP200313.AMUU09
Technical Construction File No. TPPEMJ2003113549

Certificate's Holder: XXXXXXXXXXXXXXXX
No.5 Miya Road, Yuquan Industrial Zone, Huangdong Village, Fenggang Town, Dongguan, China. 523696.

Certification ECM Mark:  Type Approved

Product: Anti Particulate Face Mask
Model[s]: AC9503, AC9003, AC8001

Verification to: Standard: EN 149:2001+A1:2009
Related to CE Directive[s]: 2016/425 (Personal Protective Equipment)


Note: The product(s) has been verified on a voluntary basis. The product(s) satisfies the requirements of the Certification Mark of ECM. If the product(s) to the above listed standards, the above Compliance Mark can be affixed on the product(s) accordingly to the ECM regulation about its release and its use. The regulation can be found at www.enfocem.com. The Certificate of Compliance can be checked for validity at www.enfocem.com (the verification does not guarantee the production of the product(s)).

Additional information: qualification about the CE marking: The holder has a 10% for the CE marking process & in place. Whereas the Manufacturer is responsible to start the CE marking Certification Procedure through the authorized notified body and to perform all the necessary activities, as required by the Directive and accepted by the notified body, before placing the CE Mark on the product(s).

Date of issue: 15 March 2020
Deputy Manager: 



Expiry date: 12 March 2025
Deputy Manager: 

Info: Certificazione Macchine Srl
Via Cof' Belletti, 243 - Loc. Castello di Senovalle - 40053 Viduggio (BO) - ITALY
☎ +39 051 4705141 ✉ +39 051 4705156 🌐 info@enfocem.it 🌐 www.enfocem.it


Ref. No. TPPEMJ2003113549

TEST REPORT

EN 149:2001+A1:2009
Respiratory protective devices - Filtering half masks to protect against particles - Requirements, testing, markings

Report Reference No.:	TPPEMJ2003113549
Compiled by (+ signature)	John Chen 
Approved by (+ signature)	Bowien 
Date of issue	2020-2-20
Number of pages	15
Testing Laboratory	Mambocert (shanghai) Technology Co., Ltd Building 17, No.912, Bibo Rd., Zhangjiang High-Tech Park, S.H.
Testing location/procedure	Mambocert (shanghai) Technology Co., Ltd Building 17, No.912, Bibo Rd., Zhangjiang High-Tech Park, S.H.
Applicant's name	XXXXXXXXXXXXXXXXXXXX
Address	No.5 Miya Road, Yuquan Industrial Zone, Huangdong Village, Fenggang Town, Dongguan, China. 523696
Test specification:	
Standard	EN 149:2001+A1:2009
Test procedure	Commission
Non-standard test method	N/A
Test Report Form No.:	EN 149:2001+A1:2009
TRF Originator	SCC
Master TRF	Dated 2006

Copyright This report is based on a blank test report prepared by Shanghai AIZEU SAFETY CO., Ltd using information obtained from the TRF originator.
Note: This report shall not be reproduced except in full, without the written approval of Shanghai AIZEU SAFETY CO., Ltd This document may be altered or revised by Shanghai AIZEU SAFETY CO., Ltd personnel only, and shall be noted in the revision section of the document.

Ref. No. TPPEMJ2003113549

Test item description	Anti Particulate Face Mask
Trade Mark	N/A
Test model/Type reference	AC9503
Manufacturer	XXXXXXXXXXXXXXXXXXXX
Factory location	No.5 Miya Road, Yuquan Industrial Zone, Huangdong Village, Fenggang Town, Dongguan, China. 523696



Summary of testing:
Ambient temperature :22°C-24°C, Humidity:50%-55%RH

Complete test was conducted on AC9503
AC9503, AC9003 are the same series products.

A representative sample of the product covered by this report has been tested and complies with the applicable requirements of this standard.