

ST3-S0110

SMARTY SAVER MINI
TRAINER SEMI-AUTOMATICO

ST7-S0110

SMARTY SAVER MINI GEO
TRAINER SEMI-AUTOMATICO

DIMENSIONI MINI
(108X119X42 - 280G)

- THE MINI TRAINER FOR SMARTY SAVER IS A TEACHING DEVICE DESIGNED TO TRAIN FIRST RESPONDERS IN THE USE OF AEDS FROM THE SMARTY SAVER RANGE, IN PARTICULAR THE SM1 (SMARTY SAVER SEMIAUTOMATIC) AND SM5 (SMARTY SAVER GEO SEMIAUTOMATIC) MODELS.



NON-SHOCKING UNIT



MOQ
3PCS

BOX CONTENT ST3-0001 / ST7-0001

SET OF UNIVERSAL REUSABLE TRAINING PADS
ACCUMULATOR WITH ITS CHARGER
REMOTE CONTROL
QUICK REFERENCE CARD
CARRYING CASE

CONTROLS AND INDICATORS

BUTTONS

- ON / OFF (KEYBOARD) (LATERAL SWITCH)
- SHOCK BUTTON: TO SIMULATE DEFIBRILLATION SHOCK DELIVERY
- ADULT MODE: TO SELECT THE PATIENT TYPE
- PAEDIATRIC MODE: TO SELECT THE PATIENT TYPE

VISUAL INDICATORS

- “AED READY FOR USE ” :1 GREEN LED
- “AED FAULT /LOW BATTERY”: 1 RED LED
- “PLACE THE DEFIBRILLATION PADS”: 2 RED LEDS
- “DO NOT TOUCH THE PATIENT: 2 RED LEDS
- “IT IS POSSIBLE TO TOUCH THE PATIENT”:1 GREEN LED
- ADULT PATIENT :1 GREEN LED
- PAEDIATRIC PATIENT: 1 GREEN LED
- ON/OFF BUTTON: 2 GREEN LEDS
- SHOCK BUTTON: 4 RED LEDS
- CPR FEEDBACK (GEO): 1 RED LED, 1 ORANGE LED, 1 YELLOW LED, 1 GREEN LED
- CPR BLUETOOTH CONNECTION:(GEO) 1 GREEN LED

SOUND INDICATORS

- MULTILINGUAL VOICES FOR INSTRUCTIONS WHILE USING THE DEVICE
- AUDIBLE WARNING AND DANGER SIGNALS

SPEAKER

- VARIABLE VOLUME THROUGH REMOTE CONTROL (POWERED WITH NO. 2 AAA BATTERIES)

DIMENSIONS & WEIGHT

- 108X119X42 MM (FOLDED HANDLE)
- 138X119X42 MM (EXTENDED HANDLE)
- 280 G (WITH BATTERY) + 23 G (ELECTRODE INTERCONNECTION CABLE)

TPOLOGY

- NO. 2 AA BATTERIES

VOLTAGE / CAPACITY

- 3V – 1AH

DURATION

- STAND-BY 5 DAYS
- CONTINUOUS OPERATION 6 HOURS